

# AUNTY




**UN**CHAINED  
LABS

## Dominate protein stability


Aunty is the fastest and highest throughput protein stability characterization tool on the planet. Period. With unrivaled fluorescence, static light scattering and dynamic light scattering skills, Aunty blasts through thermal melting and aggregation experiments – reading an entire 96-well plate of samples every minute of a thermal ramp. Sift through piles of candidates, conditions, or formulations, with as many replicates as you want. Aunty cranks out plates of results with unprecedented data resolution, so you can find your winners at breakneck speed.

- Melting temperature:  $T_m$ ,  $T_{onset}$
- Aggregation:  $T_{agg}$ ,  $T_{size}$
- Colloidal stability:  $k_D$ ,  $B_{22}$ ,  $G_{22}$
- Long term stability
- Viral vector stability




## Plate them up

The Aunty plate is the first ever 96-well quartz glass consumable for stability characterization – delivering out-of-this-world optical performance and chemical compatibility with just about any sample. Load up a plate in under a minute and snag stacks of stability experiments faster than sitting through that hour-long meeting that should've been an email. Each well needs just 8  $\mu\text{L}$  and is easily sealed tight, so your sample stays in and contaminants stay out. Prepare to fall in love with your consumable instead of struggling with it – XOXO.




## G.O.A.T. data

Aunty sets the bar for protein stability data acquisition times and quality. Its full spectrum fluorescence tracks down protein melting temperatures using either intrinsic fluorescence or reporter dyes, whichever floats your protein's boat. Add on best-in-class light scattering data to make sure you are seeing the whole picture. Aunty's dynamic and static light scattering monitor for aggregation with superior sensitivity, so you see the second it happens and how it relates to protein unfolding events. Collect oodles of data from each measurement mode separately or run them all together – Aunty is the Greatest-Of-All-Time solution for understanding the stability of your protein.




## Get eyes on aggregation

Aunty's light scattering makes sure you'll never miss aggregation in your sample. Its showstopping SLS makes for super sensitive, early detection of forming aggregates, while killer ISO compliant DLS keeps tabs on the hydrodynamic size of even the smallest number of aggregates. Aunty monitors both in parallel, so you know if aggregation is already in your prep or if it happens during a thermal ramp experiment.




## Rip through colloidal stability

Know off the bat if your protein likes its formulation or if there's aggregation on the horizon. Set up your colloidal stability experiment with Aunty and record  $k_D$ ,  $B_{22}$  and  $G_{22}$  data all in one go. Go buck wild with testing different conditions since you've got 96 wells to work with. When you're done with that experiment, send it through again at a different temperature just to see if that changes any trends.




## Speed up long-term stability

Old school long-term stability testing is clunky, clumsy and something you'd rather put off. Aunty helps you know if your protein can go the distance, while you still have time to pick a better candidate. Aunty's equipped with apps that let you monitor samples for hours or days while loaded up in the instrument. If you don't want to tie Aunty up, bring them back and forth from an incubator. Either way, your samples are safe and sound in an Aunty plate.





## Don't get left empty-handed

Viral vectors are nothing without their payload. But this valuable cargo can sneak out and ruin all the work you put into optimizing the vector in the first place. Aunty's full spectrum fluorescence is uniquely capable of tracking the escape of the genome and the stability of the capsid, so you'll know if your vector is sitting pretty or if your payload is on the lam.



## Coast through data analysis

Smash through experimental setup and get to the good stuff – actionable results. In just a few clicks, pick your app, load a sample list and let Aunty set the parameters. If you have a special experiment up your sleeve, you can take control and set all the details yourself. Now hit start and watch your results pour in live. Aunty's software lets you focus on one sample, a few, or plates of them all at once, from multiple experiments, even if they were run on different days. Overlay graphs to compare results, rapidly identify trends and draw the right conclusions in the time it takes to finish your coffee.



## Ready, set, go 21 CFR part 11

If you're working in a regulated environment, you know that controlling access to your data and keeping its integrity are must-haves. Aunty has 21CFR11 compliance software tools that help you lock your data down.

## Automate stability testing

Aunty is totally ready to integrate with lab automation. Its SBS-format 96-well plate can be loaded and moved by liquid handlers and there's an API to control getting all that data.

### 21CFR11

View your audit trail and list of experiments here.

Audit Trail Experiments Reports **Methods**

group by area Drag a field here to group by that field

	Method Name	Last Modified	Last Known Path	Integrity	Experiment Type	Created By	Updated By
7	Song and polydispersity Jun-05-2025 11553	2025Jun05 15:53:55	C:\Aunty\Methods\I	Green	Song & polydispers	Super	Super
8	Out-of-the-box stability Jun-04-2025 11505	2025Jun05 15:05:51	C:\Aunty\Methods\I	Green	Out-of-the-box stab	Supervisor2	Supervisor2
9	Out-of-the-box RD stability Jun-04-2025 115	2025Jun05 15:07:17	C:\Aunty\Methods\I	Green	Out-of-the-box stab	Supervisor2	Supervisor2
10	Isothermal Jun-04-2025 11437	2025Jun05 14:39:39	C:\Aunty\Methods\I	Green	Isothermal	Supervisor2	Supervisor2
11	Isothermal Jun-04-2025 11435	2025Jun05 14:38:39	C:\Aunty\Methods\I	Yellow	Isothermal	Supervisor2	Supervisor2
12	Isothermal Jun-04-2025 11435	2025Jun05 14:38:05	C:\Aunty\Methods\I	Green	Isothermal	Supervisor2	Supervisor2
13	Freeform Jun-09-2025 11213	2025Jun09 12:14:28	C:\Aunty\Methods\F	Green	Freeform	Super	Super
14	Freeform Jun-09-2025 11211	2025Jun09 12:12:21	C:\Aunty\Methods\F	Yellow	Freeform	Super	Super
15	Capid stability Jun-09-2025 11925	2025Jun09 09:28:14	C:\Aunty\Methods\I	Green	Capid stability	Supervisor2	Supervisor2

Export Method Details Integrity Check Back



## Specifications

Instrument	
Sample volume	8 $\mu$ L
Simultaneous samples per experiment	96
Time to read a 96-well sample plate	1 min
Temperature range	15–95 °C
Heating rate	0.1–10 °C/min
Temperature accuracy	$\pm$ 0.1 °C
Physical	43 cm W x 59 cm D x 39 cm H; 32 kg
Electrical	Auto switching power supply, 110–240 V AC, 50–60 Hz, single phase; fuse rating (F)5AH, 250 V; max. power 400 W
Laser classification	Class 1
Regulatory compliance	Software has optional 21CFR11
Fluorescence	
Parameters determined	$T_m$ , $T_{onset}$ , $T_m$ of vector payload ejection, rate of isothermal unfolding
Excitation	280 nm (UV) LED, 470 nm (blue) LED
Detection	Spectrometer at full 250–750 nm spectral range
Compatibility	Protein intrinsic fluorescence Reporter fluorescence (DSF, nucleotide dye, membrane protein dye)
Sample precision	$<1\%$ CV ( $T_m$ )
Protein sample concentration	0.025–300 mg/ml IgG (protein dependent)
AAV genome concentration	$\geq 5 \times 10^{11}$ vg/mL
Light scattering	
Parameters determined	$T_{agg}$ , $T_{size}$ , rate of isothermal aggregation Z-average diameter, size distributions, polydispersity colloidal stability: $k_D$ , $B_{22}$ , $G_{22}$
Excitation	660 nm (red) laser
Detection	Avalanche photodiode
SLS resolution	$\geq 15$ kDa change of mean MW
DLS hydrodynamic diameter range	0.3–1000 nm
DLS Size accuracy	$\pm 2\%$
Molecular weight range	192 Da–25 MDa
Minimum sample concentration	SLS: 0.05–300 mg/mL IgG (protein dependent) DLS: $\geq 0.1$ mg/mL lysozyme



**Unchained Labs**  
4747 Willow Road  
Pleasanton, CA 94588  
Phone: 1.925.587.9800  
Toll-free: 1.800.815.6384  
Email: [info@unchainedlabs.com](mailto:info@unchainedlabs.com)

© 2026 Unchained Labs. All rights reserved.  
The Unchained Labs logo, Aunty and the  
Aunty logo are trademarks and/or registered  
trademarks of Unchained Labs.