

Unchained Labs makes automation workflows a cinch, releases LEA 9!

September 19, 2017 – Pleasanton, CA – Unchained Labs, the life sciences company that’s all about getting biologists researchers the right tool for the job, released its [LEA 9 software](#) today, making [Freeslate](#) and [Junior](#) automation systems way more accessible.

LEA, short for Laboratory Execution and Analysis Software, controls Freeslate and Junior systems, executes biologists and chemistry workflow designs, integrates third-party tools, and keeps all the data in one place.

LEA 9 adds a step-by-step workflow tool called Design Creator into the mix. There’s no coding, scripting or messing around with protocols. Scientists just drag and drop their workflow steps into place and they’re good to go. Design Creator also double checks each step to make sure it’s doable, so newbie users don’t have to worry about screwing things up.

“Freeslate and Junior already help researchers wipe out major hassles from their toughest workflows,” said Taegen Clary, VP of Marketing at Unchained Labs. “Automation can be insanely complicated, and really nerve wracking for casual users. LEA 9 makes lab automation totally easy and approachable for anyone.”

About Unchained Labs

Here’s the deal. We’re all about helping biologists researchers break free from tools that just don’t cut it. Unleashing problem-tackling products that make a huge difference in the real science they do every day. That’s our mantra, our promise and we own it. We’re located in Pleasanton, CA and can be found online at www.unchainedlabs.com.

Contact:

Taegen Clary
VP of Marketing, Unchained Labs
taegen.clary@unchainedlabs.com
925.587.9806

