

Extended work surfaces

For Big Kahuna and Junior

When your current deck space won't cut it, extended work surfaces are an easy addition to any Big Kahuna or Junior to provide you that additional operating room when you need it most.

Specifications

Extended work surface (PN: SG155958, SG155958-001, F158265)



Description	Extended work surface 4 U	Extended work surface 8 U	Extended work surface for centrifuge
Unchained Labs part number	SG155958	SG155958-001	F158265
Width	270 mm	540 mm	540 mm
Depth	593 mm	593 mm	593 mm
Height	Adjustable to be coplanar with deck	Adjustable to be coplanar with deck	Coplanar with deck
Mounting location	Left or right of deck	Left or right of deck	Left or right of deck
Interface	Standard deck elements	Standard deck elements	None (centrifuge dedicated)
Element size (U) width	4 U (240 mm)	8 U (480 mm)	8 U (480 mm)
Robot accessibility	Full accessibility	Full accessibility	Full accessibility