

Rack/plate carousel

For Big Kahuna and Junior

Keep those hands off the deck. The rack/plate carousel is a functional plate and rack storage unit for use between Big Kahuna and Junior systems enabling full automation and multiple Big Kahuna systems to exchange items back and forth.



Specifications

Rack/plate carousel (PN: F157800)

Description	Specification
Physical	560 mm W x 610 mm D x 570 mm H

Model types

Unchained Labs part number	Compatible with vial racks	Number of decks	Number of total racks
F157800-001	1 mL, 2 mL	5	40
F157800-002	4 mL	4	32
F157800-003	8 mL	3	24
F157800-004	20 mL	3	24

Plates compatible with rack/plate carousel (PN: F157800)

Unchained Labs part number	Plate type
F158359	Standard vial rack
F159120	Tecan tip rack
F158972	Eppendorf PDT rack
Varies	Rainin capillary pistons
Varies	Standardized microtiter plates (including deep well)