

Vertical plate hotel

For Junior

Similar to the carousel designed for Big Kahuna, the vertical plate hotel is an element that provides an overnight stay location for your plates on a Junior. It is a vertical stack that holds microtiter plates or tip racks.

Specifications

Vertical plate hotel (PN: F161772)

Description	Specification
Physical	263 mm W x 478 mm D x 280 mm H
Height above deck	170 mm



Two different angles of the vertical plate hotel on a Junior

Maximum capacity

Capacity	Plate sizes	Maximum plate height
15 plates	<ul style="list-style-type: none"> • Shallow well plate, 370 μL (PN: V11339) • Deep well plate, 2.2 mL (PN: V139188) • 1 mL (PN: F158359) • 2 mL (PN: F157402) 	48 mm
12 plates	<ul style="list-style-type: none"> • 4 mL (PN: F157403) • 8 mL (PN: F157404) • 20 mL (PN: F157405) 	68 mm
9 plates	<ul style="list-style-type: none"> • 40 mL (PN: F157405) • 125 mL (PN: F158100) 	108 mm
9 boxes (Rainin PDT)	<ul style="list-style-type: none"> • 10 μL (PN: V148573) • 25 μL (PN: V148575) • 50 μL (PN: V148576) • 100 μL (PN: V148577) • 250 μL (PN: V148578) • 1000 μL (PN: V148579) 	111 mm
6 racks (Eppendorf PDT)	<ul style="list-style-type: none"> • 0.1 mL (PN: V15589) • 0.2 mL (PN: V16029) • 0.5 mL (PN: V16934) • 1.0 mL (PN: V15590) • 2.5 mL (PN: V15591) • 5.0 mL (PN: V15592) • 10.0 mL (PN: V122246) 	176 mm



Unchained Labs

6870 Koll Center Pkwy
 Pleasanton, CA 94566
 Phone: 1.925.587.9800
 Toll-free: 1.800.815.6384
 Email: info@unchainedlabs.com