


Wash stations

For Big Kahuna and Junior

Keep the reusable tips and needles you use on your Big Kahuna and Junior clean with a variety of wash stations, allowing you to avoid cross contamination and ensure accuracy. Each wash station comes with automated washing capabilities.

Specifications

Wash station for pH or Cavro probes with 4-channel spray (PN: F160085), wash station, multi-tip (PN: F158349) and wash station, fixed 8-tip (PN: F163541)



Description	Wash station for pH or Cavro probes with 4-channel spray	Wash station, multi-tip	Wash station, fixed 8-tip
Unchained Labs part number	F160085	F158349	F163541
Element unit (U) width	1 U	1 U	1 U
Physical	59.6 mm W x 525 mm D x 215 mm H	59.6 mm W x 525 mm D x 215 mm H	59.6 mm W x 525 mm D x 215 mm H
Vial storage positions	5 racks for 4 vials each	7	N/A
Vial maximum size for storage positions	1 and 2 mL vials	20 mL vial	N/A
Type of wash	Flush wash, spray wash and flush to drain	Flush wash or spray wash	Passive wash

Flush wash

Description	Wash station for pH or Cavro probes with 4-channel spray	Wash station, multi-tip
Tip flush wash positions	<ul style="list-style-type: none"> • Flush wells: single, 4-probes on 9 mm or 13 mm • Flushing to drain: single, 4-probes on 9mm or 13mm/ 6-probes on 9mm 	<ul style="list-style-type: none"> • Single • 3-probes on 18 mm pitch • 4-probes on 9 mm, 13 mm or 30 mm pitch • 6-probes on 9 mm, 18 mm or 20 mm pitch
Depth of flush wells	75 mm	75 mm
Maximum tip diameter	12 Ga = 2.77 mm	12 Ga = 2.77 mm

Spray wash and dry

Description	Wash station for pH or Cavro probes with 4-channel spray	Wash station, multi-tip
Pitch of tips	Single, 4-probes on 9 mm or 13 mm (manual change over)	N/A
Maximum tip diameter	3 mm	3 mm
Depth of spray well	80 mm	80 mm
Spray depth below deck	25 mm	8 mm
Normal spray flow	200 mL/min	60 mL/min
Maximum spray flow	N/A	150 mL/min
Air dry nominal flow	8 L/min	N/A



Unchained Labs

6870 Koll Center Pkwy
 Pleasanton, CA 94566
 Phone: 1.925.587.9800
 Toll-free: 1.800.815.6384
 Email: info@unchainedlabs.com

© 2018 Unchained Labs. All rights reserved. Big Kahuna and Junior are trademarks and Unchained Labs is a registered trademark of Unchained Labs. All other brands or product names mentioned are trademarks owned by their respective organizations.