

Small molecule purge box or glove box

For Big Kahuna and Junior

Air or moisture-sensitive chemistry? Inert gas purged enclosures for Big Kahuna and Junior provide a clean, low-humidity environment for moisture-sensitive work. Different styles, sizes and accessory options are available.



Specifications

Small molecule purge box or glove box

Description	Small purge box (LC Technology Solutions)	Large purge box (LC Technology Solutions)	Triple glove box (Mbraun)
Part number	602-2002	<ul style="list-style-type: none"> 602-2001 602-2003 	F160120
Compatible Systems	Junior	<ul style="list-style-type: none"> Big Kahuna (1.5 m/2.0 m) Up to two Juniors 	Big Kahuna (1.5 m/2.0 m)
Physical	309 cm W x 105.5 cm D x 222.5 cm H	312 cm W x 104.5 cm D x 222.6 cm H	590 cm W x 137 cm D x 220 cm H
Gas type	Nitrogen, argon or helium at 80 psi		
Options	<ul style="list-style-type: none"> Access door Add-on sections (121.9 cm W x 81.2 cm D) Antechambers: <ul style="list-style-type: none"> Small (15 cm D x 40 cm L) Large (39 cm D x 75 cm L) Box cooling unit Chiller feedthrough kit H₂O and O₂ Analyzer Recirculating box SpeedVac dropout Vacuum pump Ventilation system upgrade 		N/A