

# Large molecule enclosure with HEPA filter

## For Big Kahuna

Don't let air contamination hamper your experimentation. All Big Kahunas have options for HEPA enclosures providing a low bioburden arena. Each enclosure comes complete with blower, fans and filter.

### Specifications

**HEPA enclosure: 1.5 m robot (PN: F158481) and**  
**HEPA enclosure: 2.0 m robot (PN: F163108)**



Description	HEPA enclosure: 1.5 m robot	HEPA enclosure: 2.0 m robot
Physical	1,933 mm W x 1,524 mm D x 2,571 mm H	2,433 mm W x 1,524 mm D x 2,571 mm H
HEPA efficiency	99.997%	99.997%
Cleanliness class	10,000 or better in working area	10,000 or better in working area
Blower motor speed	Microprocessor controlled	Microprocessor controlled
Features	<ul style="list-style-type: none"> <li>• 2 sets sliding doors</li> <li>• Lockable/unlockable with key</li> <li>• 905 mm W x 1,170 mm H and 1,155 mm W x 1,170 mm H</li> </ul>	<ul style="list-style-type: none"> <li>• 2 sets sliding doors</li> <li>• Lockable/unlockable with key</li> <li>• 905 mm W x 1,170 mm H and 1,155 mm W x 1,170 mm H</li> </ul>