

System enclosures with ventilation

For Big Kahuna and Junior

Your science is important. Your Big Kahuna and Junior are important. Keep your systems and your science in tip-top shape with ventilated enclosures for both systems.

Specifications

1.5 m enclosure (Big Kahuna, PN: F160986), 2.0 m enclosure (Big Kahuna, PN: F160673) and ventilated enclosure (Junior, PN: AS01116)



Description	1.5 m enclosure (Big Kahuna)	2.0 m enclosure (Big Kahuna)	Ventilated enclosure (Junior)
Physical	193.3 cm W x 152.4 cm D x 218.2 cm H	243.3 cm W x 152.4 cm D x 218.2 cm H	96.8 cm W x 79.0 cm D x 117.5 cm H
Support table	300 series stainless steel	300 series stainless steel	N/A
Panels	Clear acrylic	Clear acrylic	Clear acrylic
Structure	Clear anodized aluminum, 300 series stainless steel	Clear anodized aluminum, 300 series stainless steel	Anodized aluminum
Features	<ul style="list-style-type: none"> • 2 sets sliding doors • Lockable/unlockable with key • Ventilation duct connector 	<ul style="list-style-type: none"> • 2 sets sliding doors • Lockable/unlockable with key • Ventilation duct connector 	<ul style="list-style-type: none"> • Flip panel with magnetic latch • Ventilation duct connector

Compatible accessories

Description	Carousel enclosure (Big Kahuna)	Table with side desk (Junior)
Physical	64.8 cm W x 101.6 cm D x 140 cm H	60.96 cm W x 50.8 cm D x 700–109 cm H (adjustable height)
Support table	300 series stainless steel	N/A
Panels	Clear polycarbonate	N/A
Structure	Clear anodized aluminum, 300 series stainless steel	N/A
Features	Front access through hinged panel with interlocking switch	N/A