

# Whisk


## For Big Kahuna and Junior

The whisk is an innovative arm element that helps break up aggregates and enables uniform distribution of your precursors for either the Big Kahuna or Junior.

### Specifications

#### Whisk (PN: F159963)

Description	Specification
Physical	70 mm D x 190 mm H
Maximum shaft rotational speed	120 rpm
Standard mixing container	PPR glass vial (PN: S10307-001)
Wire whisk	Titanium
Whisk drive shaft	316 stainless steel
Mounting hardware	316 stainless steel



The whisk attachment in PPR vials on the vertical shaker



Whisk attachment