

# Vertical shaker

## For Big Kahuna and Junior

Generating high intensity shaking that will fluidize solid support particles during the impregnation process, the vertical shaker is a key module to synthesize uniform catalysts with Big Kahuna or Junior.

### Specifications

#### Vertical shaker (PN: F159003)

Description	Specification
Physical	238.6 mm W x 525 mm D x 215 mm H
Plate compatibility	<ul style="list-style-type: none"> <li>• 6x8 PPR vial rack – 1x</li> <li>• 6x4 PPR vial rack – 2x</li> </ul>
Shaking amplitude	<1.5 mm
Amplitude adjustment	Manual
Frequency	50/60 Hz (same as AC frequency)

