

Crimping and decrimping tools For Big Kahuna and Junior

Automated crimping tools reproducibly work with consistent force to eliminate repetitive motion strain from the manual process. Big Kahuna and Junior use these pressure-regulated crimping and decrimping tools to maximize productivity.



Specifications

Crimping tool (PN: AS00258)

| Name | Description | Throughput (vials/hour)* | Part number |
|------------------------|----------------------------|--------------------------|-------------|
| Crimping tool 13 mm | Crimper for 13 mm seals | 290 | AS00258-13C |
| Crimping tool 20 mm | Crimper for 20 mm seals | 290 | AS00258-20C |

Decrimping tool (PN: AS00258)

| Name | Description | Throughput (vials/hour)* | Part number |
|--------------------------|------------------------------|--------------------------|-------------|
| Decrimping tool 13 mm | Decrimper for 13 mm seals | 250 | AS00258-13D |
| Decrimping tool 20 mm | Decrimper for 20 mm seals | 250 | AS00258-20D |

*Crimping/decrimping throughput will vary based on the system configuration.

Accessories for crimping and decrimping tools

Crimp seal rack

| Name | Description | Dimensions | Part number |
|-----------------------|--|--------------------------------------|-------------|
| 13 mm crimp seal rack | Crimp seal rack for 13 mm seals, 24 seals/rack | 127.76 mm W x 85.48 mm D x 36.4 mm H | AS00221 |
| 20 mm crimp seal rack | Crimp seal rack for 20 mm seals, 12 seals/rack | 127.76 mm W x 85.48 mm D x 36.4 mm H | AS00211 |

Up to 4 crimp seal racks can be vertically stacked on a deck location to maximize space.

Crimp seals

| Name | Description | Unit | Wheaton part number | Unchained Labs part number |
|-------------------|-------------------------------|---------|---------------------|----------------------------|
| 13 mm crimp seals | Aluminum center disc tear-out | 1000/cs | 224192-01 | V17755 |
| 20 mm crimp seals | Aluminum center disc tear-out | 1000/cs | 224193-01 | V17756 |

Other crimp seals and materials may be compatible.



Unchained Labs

6870 Koll Center Pkwy
Pleasanton, CA 94566
Phone: 1.925.587.9800
Toll-free: 1.800.815.6384
Email: info@unchainedlabs.com

© 2018 Unchained Labs. All rights reserved. Big Kahuna and Junior are trademarks and Unchained Labs is a registered trademark of Unchained Labs. All other brands or product names mentioned are trademarks owned by their respective organizations.