

# Cooled storage bay

## For Big Kahuna and Junior

We know keeping a consistent temperature is key to your screening success. Big Kahuna and Junior have a cooled storage bay module that includes three independently cooled zones keeping your formulations nice and chilly.

### Specifications

Cooled storage bay (PN: F158763)

Description	Specification
Element size (U) width	2 U (120 mm)
Plate size	Microtiter footprint (127.76 mm W x 84.48 mm D x 25.4 mm H)
Maximum plate height	122 mm
Number of fully independent cooled zones	3
Vial operating temperature*	Ambient - 4 °C
Vial temperature accuracy	±2 °C
Vial temperature uniformity within each cooled zone	2 °C
Cooling ramp rate	0.7 °C/min

\*Vial operating temperature is vial- and solution-dependent.

